

SIMPATI version	Desktop OS								Server OS				
	Win10	Win8.1	Win8	Win7 <sup>1)</sup>	Vista	WinXP	Win200	Win98	2016 <sup>3)</sup>	2012	2008 <sup>1)</sup>	2003	NT4
2016 <sup>2)</sup>	●	●	●	●					●	●	●		
4.50SP2		●	●	●	●					●	●		
4.50SP1		●	●	●	●						●		
4.50		●	●	●	●						●		
4.06		● <sup>4)</sup>	● <sup>4)</sup>	●	●	●					●	●	
4.04				●	●	●					●	●	
4.02						●						●	
4.0						●						●	
3.0						●	●					●	
2.xx						●	●	●				●	●
1.xx								●					●

<sup>1)</sup> 32Bit or 64Bit versions

<sup>2)</sup> current version (03/2018)

<sup>3)</sup> the SIMPATI PC may not simultaneously be configured as domain controller

<sup>4)</sup> only with installed patch 1+2 for 4.06